

Pwo DNA Polymerase

Specification

Catalog No.	1654 101 Thermostable DNA polymerase
Origin	Pyrococcus woesei, recombinant in E. coli RM82
Appearance	clear, colorless solution in Tris-HCl, 20 mmol/l; KCl, 100 mmol/l; EDTA, 0.1 mmol/l; DTE, 1 mmol/l; Nonidet P40, 0.5% (v/v); Tween 20, 0.5% (v/v); glycerol, 50% (v/v); pH approx. 7.5 at 20°C
Volume activity	$\geq 5 \times 10^3$ U/ml
Endonucleases (λ DNA) (MWM II) Nicking activity (pBR322)	not detectable (≤ 30 U/4 h/72°C) not detectable (≤ 30 U/4 h/72°C) not detectable (≤ 30 U/4 h/72°C)
Function testing in PCR (λ DNA, 0.5 kb) (human genomic DNA, 1.1 kb)	corresponds corresponds
Stability	at -20°C 18 months

Description

Thermostable Pwo DNA Polymerase was originally isolated from the hyperthermophilic archaebacterium *Pyrococcus woesei*. It is now supplied as a recombinant enzyme from *E. coli*. Pwo DNA Polymerase has a molecular weight of about 90 kD. It is a 5' – 3' DNA polymerase and has a 3' – 5' exonuclease (proofreading) activity. The enzyme has no detectable 5' – 3' exonuclease activity.



Application

Because of its proofreading activity, Pwo DNA Polymerase amplifies DNA with tenfold greater accuracy than other common thermostable polymerases. It is especially useful for such demanding applications as

- Routine high fidelity amplification of up to 3 kb human genomic DNA targets.
- Amplification of genomic products for cloning.
- Study of allelic polymorphisms.
- Characterization of rare mutations.
- Characterization of cell populations.
- Labeling of PCR products with modified nucleotides (DIG-dUTP, biotin-dUTP, fluorescein-dUTP).

Pwo DNA Polymerase cannot amplify targets longer than approx. 3 kb. For high fidelity amplification of longer templates, see the product profiles of the Expand PCR Systems.

Key advantages

- Maximum PCR fidelity, because the proofreading activity of Pwo DNA Polymerase excises misincorporated nucleotides.
- High PCR yield, because the polymerase is very resistant to the high denaturing temperatures used in PCR.
- Simplifies cloning of PCR products because it generates blunt-ended products.

Properties

Thermal stability	Pwo DNA Polymerase is stable for > 2 h at 100°C.
Divalent ion requirement	Mg ²⁺ (standard concentration, 2 mmol/l). The polymerase prefers MgSO ₄ to MgCl ₂ . The enzyme is particularly sensitive to Mg ²⁺ concentration when the target is > 2 kb long.
dNTP concentration	approx. 200 µmol/l for each dNTP. Caution: Lower deoxynucleotide concentrations can activate the 3' – 5' exonuclease proofreading activity, which might degrade primers and products.
Effect of additives	Usually, additives will not improve the performance of Pwo DNA Polymerase in PCR. Nevertheless, additions of up to 100 µg/µl bovine serum albumin (BSA), Triton X-100 (0.1%), pyrophosphatase, DMSO, or betaine may improve efficiency in some cases.